

AquaRNA Instruction Manual

General Information

Description

AquaRNA™ is a multifunctional aqueous solution-based reagent for DNA, RNA, and protein extraction. This single solution will lyse the cells, inactivate degradative enzymes, and extract DNA, RNA, and proteins. DNA and RNA are recovered from the cell lysate by isopropanol precipitation, while proteins remain soluble in the AquaRNA-isopropanol solution and can be recovered by acetone precipitation. AquaRNA enables concurrent isolation of DNA, RNA, and proteins from the same specimen without using different DNA, RNA, and protein extraction kits.

Specification

Product Name	AquaRNA™ Kit
Product #	5001, 5030
Size	5001: 1 ml; 5030: 30 ml
Kit Contents	5001: 1 ml AquaRNA Solution, User Manual 5030: 30 ml AquaRNA Solution, User Manual
MSDS	Available at www.aquaplasamid.com
Storage	Store tightly capped at 22 °C. Vortex to mix well before dispensing.
Note	In addition to AquaRNA, please order ProSink (# 9030) for protein removal when extracting RNA from blood, liver, and spleen; and ProMelt (# 1115) for protein solubilization when extracting proteins.

Terms & Condition

Product Usage: For In Vitro Laboratory Research Use Only. NOT to be administered to humans or used for medical diagnosis.

Limited Product Warranty: We offer a LIMITED PRODUCT WARRANTY to our customers. This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by MultiTarget Pharmaceuticals. We shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Product Warning: Contains guanidine thiocyanate. Harmful if swallowed. Causes irritation to skin, eyes, and respiratory tract. Do not mix with Bleach.

Patents, Trademarks & Copyrights

AquaRNA™ is a trademark of MultiTarget Pharmaceuticals, LLC. ©2017 Multitarget Pharmaceuticals, LLC. All rights reserved.

AquaRNA Cell Protocol

This protocol uses 100 µl of AquaRNA to extract DNA, RNA, and proteins from 0.5-2 million cultured mammalian cells or 0.25 ml log-phase microbial culture.

1. Harvest the cells

For eukaryotic cells: Centrifuge at 3,000 xg for 5 min to pellet 0.5-2 million cultured cells in a 1.5-ml microfuge tube, and aspirate to remove the medium.

For microbial cells: Centrifuge 0.25 ml log-phase culture at 12,000 xg for 1 min to pellet the cells. Aspirate to remove the medium and suspend the cells in 100 µl of 1 mg/ml lysozyme (not included, use lyticase or equivalent for yeast cells) in TE buffer containing 300 mM KCl (*pH 8-9, lysozyme will not be as effective at pH <8*) and incubate on ice for >15 minutes, and vortex occasionally. Centrifuge to pellet the cells and aspirate to remove the supernatant.

2. Extract the DNA/RNA

Add 100 µl of AquaRNA to the cell pellet. Vortex to mix well (*for Gram+ bacteria, use a bead beater to assist cell lysis*). Incubate at room temperature for 5min. Vortex and centrifuge at 12,000 xg for 5 min to pellet the cell debris.

3. Pellet the DNA/RNA

Transfer the clear lysate (80 µl) to a new 0.5-ml microfuge tube. (A) To pellet total DNA/RNA: Add 0.7 vol (56 µl) of isopropanol to the lysate and vortex to mix. Centrifuge at 12,000 xg for 5 min to pellet the DNA/RNA. Decant to discard the supernatant (*Proteins remain in the isopropanol supernatant and can be recovered by precipitation in 4 vol of acetone or by dialysis*). (B) To pellet RNA and DNA differentially: Add 0.25 vol (20 µl) of isopropanol to the lysate, vortex and centrifuge at 12,000 xg for 5 min to pellet large RNA. Transfer the supernatant containing DNA and small RNA to a new tube. Add an additional 40 µl of isopropanol, vortex and centrifuge for 5 min to pellet the DNA and small RNA. To rinse the DNA or RNA pellet: Gently fill the tube with 70% ethanol from a squirt bottle and then decant to discard the ethanol solution. Be sure to rinse the entire interior of the tube, including the inside of the cap. Repeat the ethanol rinse once (use centrifuge if the pellet is loosening). Tap the tube on a paper towel to remove residual ethanol, place the tube upside down on the paper towel to air dry the DNA/RNA pellet for 5-10 min. Add 100 µl of nuclease-free water to the pellet, vortex and/or pipette to solubilize the DNA/RNA pellet. Incubate at 22 °C for 5 min, centrifuge at 12,000 xg for 5 min to pellet any insoluble material, and transfer the DNA/RNA solution to a new tube.

AquaRNA Tissue Protocol

This protocol uses 0.5 ml of AquaRNA to extract DNA, RNA, and proteins from 25 mg of animal tissues or 50 mg of plant tissues.

1. Homogenize the tissue

Homogenize the animal tissue (~25 mg) or plant tissue (~50 mg) in 0.5 ml AquaRNA with a pestle-and-tube homogenizer (or a multi-channel bead beater). Move the pestle up and down while vortexing at top speed to fully homogenize the tissue. Add 2-3 drops (20-30 μ l) of isopropanol to the homogenate to reduce foaming, push the pestle to the bottom of the homogenizer and pour the homogenate into a 1.5-ml microfuge tube. Centrifuge at 12,000 xg for 5 min to pellet the tissue debris (*For plant samples, add ~0.25 vol of AquaRemove (#1208, order separately and dilute it with an equal volume of isopropanol before use), e.g., add 100 μ l AquaRemove to 300 μ l of the crude lysate, invert to mix well, incubate at -20 °C for 5 min and centrifuge again to pellet the insoluble*).

2. Recover the DNA/RNA

Transfer the clear lysate (~400 μ l) to a new 1.5-ml microfuge tube. Add 0.7 vol (~280 μ l) of isopropanol and vortex to mix. Centrifuge at 12,000 xg for 5 min to pellet the DNA/RNA. Decant to discard the supernatant (*Proteins remain in the isopropanol supernatant and can be recovered by precipitation in 4 vol of acetone as described in Step 3. Recover the proteins or by dialysis*). Gently fill the tube with 70% ethanol from a squirt bottle and then decant to discard the ethanol solution. Be sure to rinse the entire interior of the tube, including the inside of the cap. Repeat the ethanol rinse once. Tap the tube on a paper towel to remove residual ethanol, place the tube upside down on the paper towel to air dry the DNA/RNA pellet for 5-10 min. Add 400 μ l of nuclease-free water to the pellet, vortex and/or pipette to solubilize the DNA/RNA pellet. Incubate at 22 °C for 5 min, vortex and centrifuge at 12,000 xg for 5 min to pellet any insoluble material, and transfer the DNA/RNA solution to a new tube and store at -20 °C.

3. Recover the Proteins

Transfer 300 μ l protein-containing isopropanol supernatant to a 2-ml microfuge tube. Add 4 vol (1.2 ml) acetone, vortex and centrifuge at 12,000 xg for 5 min to pellet the proteins. Decant to discard the supernatant, tap the tube on a paper towel to remove residual acetone. Immediately add 400 μ l of ProMelt (#1115, order separately) to the wet protein pellet, pipette and vortex to suspend the pellet. Incubate at 22 °C for 15 min to solubilize the proteins. Vortex and centrifuge at 12,000 xg for 5 min to pellet any insoluble material. Transfer the protein solution to a new microfuge tube and use it directly in SDS-PAGE or store it at 4 or -20 °C (*Some SDS may precipitate out at these temperatures, it may be re-solubilized by incubating at 65 °C for 5-10 min*).

Frequently Asked Questions

Please read through these questions carefully. The answers provide additional helpful tips and useful information for the successful use of AquaRNA.

1. How should I store the AquaRNA solution?

It may be stored at 22 °C for 12 months. If AquaRNA becomes precipitated when exposed to low temperatures, you may incubate it at 37 °C for 15-20 min to resolubilize the reagent.

2. My RNA was degraded, where was the RNase coming from?

To troubleshoot RNase contamination, you may set up a DNase I digestion in 1x DNase buffer. Before adding DNase I, divide the sample into two aliquots and add DNase I to one of them. If RNA degradation is seen only in the DNase I treated sample, your DNase I may be contaminated. If RNA is degraded without adding DNase I, your RNA sample may be contaminated. A good habit to prevent RNase contamination is to ensure that your gloves or fingers do not touch the inside of the lid and the mouth of the tube when opening and closing the tube containing RNA solution.

3. How should I remove the genomic DNA from my DNA/RNA preparation?

You may add 0.2 units of DNase I to 10-20 µl of DNA/RNA solution in 1x DNase buffer, and incubate at 22-37 °C for 20-30 min. Then run the digested sample in a 0.8% native agarose gel to confirm that the DNA digestion is complete and the RNA bands are discrete. To inactivate the DNase I, use Ambion's DNase removal reagent or heat-inactivate the DNase I at 65 °C for 15 min.

4. Why did my DNA/RNA solution show a strong absorption below A260?

It is likely due to trace amount of guanidine salt contamination. If it interferes with your downstream applications, you may further purify the extracted DNA/RNA with a silica spin column (e.g., a plasmid miniprep column). Simply add an equal volume of 4 M GuHCl and 1M NaOAc (pH unadjusted, ~7.0) to your DNA/RNA solution (may contain the insoluble pellet) and load it into the spin column, centrifuge to allow DNA/RNA binding to the silica membrane, wash the column with 0.6 ml 75% EtOH, and elute the DNA/RNA in 50 µl of deionized water or TE buffer.

5. Can I do RT-PCR without removing the contaminating genomic DNA?

Complete DNA removal may be difficult to achieve and unnecessary if you use intron-spanning primers for the PCR amplification. You may also design and use a 5' tailed RT primer to make the cDNA and then use a pair of PCR primers with one of them complementary to the unique tailed region of the RT primer to amplify the cDNA [Hurteau and Spivack. mRNA-specific reverse transcription-polymerase chain reaction from human tissue extracts. *Anal Biochem.* 2002 Aug 15;307(2):304-15; and Chen, et al. Real-time quantification of microRNAs by stem-loop RT-PCR. *Nucleic Acids Research* 2005 33(20):e179], especially when intron-spanning is unavailable. In any case, you should always include a no-RT control in your amplification to confirm that your primers do not amplify the contaminating genomic DNA.